lopital®
If caregivers have a choice

Tango XL

Extra space, comfort and stability for higher care demands

5100 5750

Tango XL shower toilet chair

The Tango XL shower toilet chair is specially developed for clients who need extra space, comfort, and secure support during daily care routines. With its reinforced construction and generous seating dimensions, the chair provides a safe and comfortable solution for users with a higher body weight—while maintaining smooth manoeuvrability for caregivers.

Designed for use in professional healthcare environments, the Tango XL combines practical functionality with ergonomic comfort. Foldable armrests and swing-away footrests support easy transfers, while the hygienic removable seat and smooth-rolling castors with brakes ensure safe and efficient handling during showering and toilet use.

A reliable choice for care facilities that aim to deliver inclusive care with maximum comfort, durability, and hygiene.

Highly Stable Design. The robust frame, low center of gravity and braked castors ensure maximum stability during transfers and care moments.

Ergonomic Use. Equipped with smart handle and four braked wheels to simplify maneuovering, ensuring high stability on various surfaces and minimizing the strain on the caregiver's joints.

Reliable & Long-Lasting. Built for daily care and backed by a 60-month warranty, offering a dependable performance.

Antibacterial & Corrosion-Resistant. The antibacterial coating with silver-ion supports hygiene and infection control, while the corrosion-resistant frame with smooth surfaces ensures durability and easy cleaning.


Discover our tools in action through the Interactive Product videos on our website: www.lopital.com



Technical specifications	
Max. load	250 kg / 551.2 lbs
Product weight	27 kg / 59.2 lbs
Seat width	55 cm / 21.7 inch
Seat depth	49 cm / 19.3 inch
Seat height	53 cm / 20.9 inch



Suitable for people with very slightly to slightly reduced mobility



Ergonomical seat for toilet use and bedpan.



Comfort-optimized handle for smooth and light operation

Foldable armrests for improved accessibility

Ergonomic seating comfort with removable seat sections featuring a snap-lock system for easy cleaning


Guide rails for bedpan or toilet bucket

Bumper to prevent damage to the frame

Foldable, height-adjustable, swing-away footrests for smooth and safe transfers

Optional

6950 1104
Plastic bedpan
and lid.



Options:
- 6840 5500 Adjustable head/neck support
- 6950 1104 Bedpan with lid plastic
- 6975 8044 Seatbelt 2-point
- 6975 8045 Seatbelt 4-point
- 6960 8026 Comfort leg support

Technical specifications	
Total lenth base frame	106 cm / 41.7 inch
Total width base frame	65 cm / 25.6 inch
Turning radius	106 cm / 41.7 inch

For more information also visit our website: www.lopital.com
Lopital Nederland B.V., Laarakkerweg 9, 5061 JR Oisterwijk
Tel +31 (0)13 52.39.300, E-mail info@lopital.nl, www.lopital.nl

Green Deal

ISO CERTIFIED 9001 - 13485 - 14001

MADE IN THE NETHERLANDS

Scan for information

LFT5750-2026-01 EN

All of the aids produced by Lopital comply with the applicable standards for Class I medical products (MDR (EU) 2017/745).

Copyright ©2026 Lopital. All rights reserved.